Trace analysis has never been that easy!

EASE OF USE

TRUSTFUL QUALITY

The trace SN cube is optimized to significantly

easy accessible system components as well as a

ensures reliable leak-tightness of the instrument

meet the highest quality standards and reliability. They are certified and validated in accordance

SAMPLE	SULFUR [mg/kg]	RSD [%]	NITROGEN [mg/kg]	RSD [%]
DIESEL	348	0.74	166	0.65
MINERAL OIL	0.664	0.96	6.66	1.88
TECHNICAL OIL	3.84	2.23	0.06	4.53
BENZENE	0.225	0.48	0.22	4.65
KETONE	183	1.04	7.56	1.85
AQUEOUS PROTEIN SOLUTION	26.9	1.34	363	2.42
BIOMASS	27.8	1.83	385	1.31

IDEAL SOLUTION FOR

- R&D petrochemical laboratories
- Quality control laboratories



High sensitivity

Outstanding sensitivity thanks to high performance, state-of-theart technology.



High data quality

Outstanding precision and accuracy through high performance combustion. Matrixindependent results. Longterm stability of calibration.



SAMPLE TYPES ANALYZED

• LPG

Great flexibility

Wide range of optional conversion kits available for special applications. Upgradeable at any time.



Extreme durability

Outstanding robustness and longevity thanks to state-of-the-art technology. 10 year warranty on furnace and TCD cell.



Elementargroup - your partner for elemental analysis

Elementargroup is the world leader in high performance analysis of organic elements. Continuous innovation, creative solutions and comprehensive support form the foundation of the Elementar and Isoprime brands ensuring our products continue to advance science across agronomy, chemical, environmental, energy, materials and forensics markets in more than 80 countries.

Elementar Analysensysteme GmbH | Donaustraße 7 · 63452 Hanau (Germany) Phone: +49 (0) 6181 9100-0 | info@elementar.de | www.elementar.de

elementar EXCELLENCE IN ELEMENT

elementargroup.



ppb analysis at its best



trace SN 👘 cube





Great flexibility



trace SN 🍘 cube



trace SN Cube Outstanding

sample flexibility

KEY FEATURES

- Industry-leading performance

- for sulfur

- For liquids, LPG, gas and
- solid samples
- Outstanding robustness
- Up to 60 position autosampler

High sensitive detection

High sensitive chemiluminescence and UV fluorescence detectors ensure utmost sensitivity: 6 ppb fur sulfur and 15 ppb for nitrogen. This ultra-trace detection power gives rise to uncompromised analysis of sulfur and nitrogen even in samples of lowest concentrations. This ensures that the trace SN cube greatly exceeds all requirements of international standards and gives confidence to fulfill expected tighter standard requirements in future.

Outstanding sample flexibility

Along with liquid injections the trace SN cube can be equipped with gas, LPG or solids module for a variety of matrices. A temperature controllable autosampler is also available which allows the analysis of highly volatile samples. In addition, the instrument can be easily retrofitted to measure liquid samples.

The trace SN cube is the world's most sensitive combustion analyzer for sulfur and nitrogen in diesel, gasoline and related samples. The catalyst-supported digestion technology ensures complete digestion of all kinds of samples without soot formation or insufficient recovery. The sample is always injected directly into the combustion tube, there is no

need for time-consuming boat inlet feeding. Thus, the analysis time is just 3-5 minutes, depending on sample type. The trace SN cube works as two analyzers in one box. The sample is injected twice in different combustion tubes which allows the specific optimization of the oxidation process for each element.



③ Reliable determination of 100 µg / kg sulfur in a standard thanks to an outstanding low baseline noise and drift.

Unique N Excess module

The known problem of interference in UV fluorescence detection by NO (higher N concentrations simulate too high sulfur contents) is solved by means of Elementar's unique N Excess module. SO, is collected at an adsorption column and released after the decay of the signal caused by NO. This technique allows the precise determination of 0.1 mg / kg S in the presence of more than 100 mg / kg N.

True two channel operation

The trace SN cube works with two independent furnaces and combustion tubes. Thus, it is possible to adjust all parameters which may influence the results, such as temperature, injection volume, injection speed etc. of both channels completely independently. This results in an outstanding sample flexibility at always highest precision.

OIL ANALYSIS BY HIGH TEMPERATURE COMBUSTION



High temperature combustion for the still the matter of choice if sensitivity really means of UV fluorescence detection for sulfur and chemiluminescence for nitrogen determination detection limits of as low as 10 ppb can be achieved without problems, outperforming all competing techniques.



CATALYST SUPPORTED COMBUSTION

Almost any samples can be ignited without any soot production at injection volumes up to 100 µl. This offers the customer a flexibility in sample handling which is